

Dolphin Imaging Runs on a Mac!

Dolphin Imaging workstations and standalone stations can now run on recent model Macs, using either free Apple Boot Camp software, or Parallels software from Parallels Inc. OS X VNC remote control software is also required. The Dolphin Development Team still recommends using a PC as your server

- **Boot Camp.** Boot Camp is free software from Apple that allows dual booting of Mac OS X and Windows XP SP2. Only one operating system can be run at a given time. Visit <http://www.apple.com/macosx/bootcamp/> for further details.
- **Parallels Desktop.** Parallels Desktop is Virtual Machine Software from Parallels, Inc. It runs a virtualized Windows operating system simultaneously with your native Macintosh operating system. Visit <http://www.parallels.com/en/products/desktop/> for further details.
- **OS X VNC Software.** OS X VNC is remote control software that lets our tech support staff access and work inside your computer when you need technical assistance. Visit <http://www.redstonesoftware.com/products/vine/server/vineosx/index.html> to learn more about this technology.

Hardware and software requirements for running Dolphin on a Mac:

- Intel-based Macintosh Computer
- OS X 10.4.6 software or better
- Parallels Desktop and/or Apple Boot Camp
- Licensed retail copy of Windows XP Pro with SP2
- OS X VNC Remote Control Software
- Broadband Internet connection
- 2-button mouse

Intel-based Macintosh Computers:

- The Mac Mini
- The iMac
- The MacBook
- The MacBook Pro
- The MacPro

NOTE: High Resolution Dolphin 3D rendering is NOT available for the following models:

- Mac Mini
- 17-inch iMac
- standard 20-inch iMac
- MacBook.